Periform[®]+

Instruction Manual – GB

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Manual de Instrucciones – ES

Handleiding – NL

Bruksanvisning – SW

Bruerveiledning – NOR

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inner strength, outer confidence

Indications for Use:

For use as an intra-vaginal probe in the treatment of female 'Stress', 'Urge' or 'Mixed' urinary incontinence with muscle stimulation and EMG biofeedback equipment approved to EN60601-01.

**FOR SINGLE PATIENT USE – May be re-used by the same patient.
Use as directed by your medical adviser**

Contraindications & Precautions:

DO NOT USE:

With electrical stimulation if fitted with a cardiac pacemaker • If pregnant
• If the package seal is not intact • If an infection or tissue damage or recent surgical scar is present • During menstruation • With electrical stimulation in patients with a history of carcinoma at the site of stimulation • If an irritation or discomfort occurs and seek medical guidance.

CAUTION IS ADVISED:

When using electrical stimulation on patients with altered or reduced tissue sensation • During electrical stimulation if the user has a history of epilepsy - seek informed medical guidance • When using the PERIFORM⁺ with equipment capable of delivering outputs in excess of 2mA/cm². - if in doubt contact the equipment manufacturer or your medical adviser.


Cleaning:

- Wash PERIFORM⁺ in warm soapy water.
- Rinse in clean water and dry.
- Repeat after each use and dry thoroughly before storing.
- DO NOT use boiling water.

Disposal:

Clean thoroughly and dispose of according to local regulations.

Figure 1

Instructions for Use:

IMPORTANT:

- Before first and subsequent use wash and dry as directed.
- Apply light coating of KY Gel or similar lubricating gel to tip and metal surfaces of probe to aid insertion and provide good electrode conductivity.
- If using with Contraction Indicator connect together as shown in Figure 1.
- Insert into vagina in 'East/West' position, i.e. electrode surface facing towards hips; the external 'flange' should sit comfortably between the labia.
- Ensure equipment is turned off before connecting to probe cables.
- Use as directed by your medical adviser. After use, turn off equipment, disconnect cables. Do not pull on cables to remove probe.
- Wash, dry and store as instructed.


Contraction Indicator:

When fitted to the PERIFORM⁺ the Contraction Indicator provides valuable evidence of the performance of the pelvic floor muscles. It may be used on its own or when using stimulation and biofeedback equipment. When pelvic floor muscles contract correctly the indicator will move downwards ().

Moving upwards indicates incorrect contraction (X) as shown in Figure 2.

Functions and benefits of the Indicator are:

- Quantifies the strength and duration of voluntary and involuntary (stimulated) contractions. N.B. Voluntary contractions can be used to reinforce stimulated contractions.
- Indicates need to use stimulation to strengthen a weak pelvic floor – i.e. little or no indicator movement.
- Teaches bracing techniques to help avoid leakage when coughing / sneezing / straining etc.
- Demonstrates the success of a pelvic floor exercise programme

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